



Material Inovys Evo			
No.	Plan	Article	Quantity
1	Y011362	1658154 (40)	1
	Y011932	1663554 (45)	1
	Y011936	1663654 (50)	1
	YP06990	2743654 (35)	1
2	Y011922	1658354 (R)	1
3	Y011422	1658254 (L)	1
4	M8 DIN126	1996214	10
5	M8 DIN 985	2994009	8
6	M12 (PA6.6)	1916509	4
7	M8 (PA6.6)	2994415	4
8	M8x60	2995491	2
9	M8x55	1919452	2

INOVYS PIVOT FRAME				Plan No.
Sheet 1/2				E000383 / H
Scale: 2:5	Drafted	Verified	Confirmed	
Date	05/07/2007	17/05/2024	17/05/2024	
Init.	SD	PM	AS	

E000383 / H

Material Inovys II			
No.	Plan	Article	Quantity
1	Y011362	1658154 (40)	1
	Y011932	1663554 (45)	1
	Y011936	1663654 (50)	1
	YP06990	2743654 (35)	1
2	Y011922	1658354 (R)	1
3	Y011422	1658254 (L)	1
4	M8 DIN126	1001554	10
5	M8 DIN 985	1001682	8
6	M12 (PA6.6)	1916509	4
7	M8 (PA6.6)	1003199	4
8	M8x60	1900809	2
9	M8x55	1995854	2

E000383 / H

INOVYS PIVOT FRAME				Plan No.
Sheet 2/2				E000383 / H
Scale: 2:5	Drafted	Verified	Confirmed	
Date	05/07/2007	17/05/2024	17/05/2024	
Init.	SD	PM	AS	