



PART NAME	PART NUMBER
1 FRAME LEFT	1923037
1 FRAME RIGTH	1923038
2 BACK FRAME	1923039
3 ARMREST L	1923041
3 ARMREST R	1923040
4 ARMREST SUPPORT L	1923042
4 ARMREST SUPPORT R	1923043
5 FOOTREST ROD 2 PC	1923044
6 FOOTPLATE	1923045
7 FRONT FRAME	1923046
8 LITH BATTERY + HOLDER	1923047
9 REMOTE	1919239
10 REMOTE BRACKET	1923049
11 LOCK HANDLE	1923060
12 FRONT FORK	1923050
13 FRONT WHEEL	1923051
14 MOTOR LEFT	1923062
14 MOTOR RIGHT	1923063
15 REAR WHEEL	1923052
16 ANTI TIP WHEEL	1923061
17 CONTROLLER	1919238
17 CONTROLLER SEALING KIT	1919240
18 SEAT FRAME	1923053
19 SEAT PLATE	1923054
20 SEAT CUSHION	1923055
21 UPHOLSTERY BACK	1923056
22 ARMREST FOLDING ELEMENT L	1923057
22 ARMREST FOLDING ELEMENT R	1923058
23 REAR CONNECTION ROD	1923079
NOT DEPICTED BATTERY CHARGER	1923059