



Bill of Material			
No.	Plan	Article	Quantity
1	P00056945	2676154	2
2	M8x16 DIN 912	1915411	2
3	M8	1915947	2
4	P00065028	2679000	2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP

ATTACHING THE BELTS TO THE INOVYS 2 SEAT				---
				Plan No. E00065025/D.1
Scale: 1:3	Drafted	Verified	Confirmed	
Date	2023-07-22	11/12/2023	11/12/2023	
Initials	AS	PM	AS	