



| | | | | |
|--|---|--|--|--|
| SIDE FRAME ECLIPS X4/KIDS 30° ZIJFRAME ECLIPS X4/KIDS 30° | | | | Plan No. / Plannr. E000470/F |
| Scale / Schaal: 1:5 | Drafted / Getekend 29/07/2008 | Verified / Nagezien 21/10/2015 | Confirmed / Bevestigd 21/10/2015 | |
| Date / Datum 29/07/2008 | 21/10/2015 | 21/10/2015 | | |
| Init. / Init. SD | TT | PJ | | |

| Material / Materiaal | | | | | | | | | | | |
|----------------------|---------|-------------------|-----------------|-----------|---------|-------------------|-----------------|-----------|---------------------|-------------------|-----------------|
| No. / Nr. | Plan | Article / Artikel | Number / Aantal | No. / Nr. | Plan | Article / Artikel | Number / Aantal | No. / Nr. | Plan | Article / Artikel | Number / Aantal |
| 1 | X000106 | G003218 | 1 | 27 | Y004178 | 1504200 (Le/Li) | 1 | 50 | --- | 1004550 | 1 |
| 2 | X000858 | 1006467 | 1 | | Y004326 | 1504900 (Ri/Re) | 1 | 51 | --- | 1004782 | 1 |
| 3 | X000862 | 1006463 | 1 | 28 | --- | 18265.. | 1 | 52 | M6x35 | 1006057 | 2 |
| 4 | X000989 | 1004539 | 3 | 29 | Y005293 | 15502.. | 1 | 53 | --- | 1006110 | 1 |
| 5 | X001016 | 1001682 | 1 | 30 | Y005322 | 15403.. | 1 | 54 | --- | 1006111 | 4 |
| 6 | X001119 | 1001554 | 2 | 31 | Y005466 | 15470.. | 1 | 55 | --- | 1006112 | 1 |
| 7 | X001674 | 1990601 | 2 | 32 | Y005806 | 15340.. (Le/Li) | 1 | 56 | --- | 1006850 | 3 |
| 8 | X001842 | G053400 | 2 | | Y005956 | 15341.. (Ri/Re) | 1 | 57 | --- | 1006999 | 3 |
| 9 | X002340 | 1504500 | 2 | 33 | Y006196 | 1595082 | 1 | 58 | --- | 1007000 | 2 |
| 10 | X002135 | 1992405 | 1 | 34 | Y006648 | 1595182 | 1 | 59 | --- | 1008459 | 1 |
| 11 | X002205 | 1489300 | 2 | 35 | Y013003 | 16806.. | 1 | 60 | --- | 1008492 | 16 |
| 12 | X002276 | 1497400 | 1 | 36 | Y013029 | 16969.. (Le/Li) | 1 | 61 | --- | 1009588 | 1 |
| 13 | X002277 | 1992933 | 1 | | Y013035 | 16970.. (Ri/Re) | 1 | 62 | --- | 1009902 | 1 |
| 14 | X002495 | 1993677 | 1 | 37 | Y013032 | 16968.. (Le/Li) | 1 | 63 | --- | 1990648 | 3 |
| 15 | X002539 | 1993795 | 1 | | Y013033 | 16971.. (Ri/Re) | 1 | 64 | --- | 1990712 | 2 |
| 16 | X002540 | 1993796 | 1 | 38 | --- | 1001534 | 3 | 65 | M6x45 DIN912/8.8 | 1003109 | 2 |
| 17 | X002568 | 15339.. | 1 | 39 | --- | 1001561 | 2 | 66 | --- | 1992559 | 2 |
| 18 | X002636 | 1994438 | 1 | 40 | --- | 1901017 | 1 | 67 | --- | F002907 | 1 |
| 19 | X002670 | 1547400 | 1 | 41 | --- | 1001564 | 12 | 68 | M16 type S DIN17222 | 1911653 | 1 |
| 20 | X002749 | 1553400 | 1 | 42 | --- | 1001580 | 2 | 69 | --- | 1994605 | 1 |
| 21 | X003303 | 1596382 (Le/Li) | 1 | 43 | --- | 1001589 | 1 | 70 | X003481 | 1996587 | 2 |
| | X003302 | 1596282 (Ri/Re) | 1 | 44 | / | / | / | 71 | X003480 | 1996670 | 1 |
| 22 | X003384 | 1997149 | 1 | 45 | --- | 1001610 | 1 | 72 | --- | F002824 | 2 |
| 23 | X003852 | 16700 | 1 | 46 | --- | 1001670 | 2 | 73 | M6x35 DIN912/8.8 | 1001496 | 4 |
| 24 | X004006 | 1901672 | 1 | 47 | --- | 1001680 | 2 | | | | |
| 25 | Y004176 | 1504300 | 1 | 48 | --- | 1001684 | 2 | | | | |
| 26 | Y004177 | 1504400 | 1 | 49 | --- | 1002994 | 2 | | | | |

| | | | | |
|-------------------------------|--------------------|---------------------|-----------------------|---|
| SIDE FRAME ECLIPS X4/KIDS 30° | | | | Plan No. / Plannr. |
| ZIJFRAME ECLIPS X4/KIDS 30° | | | | E000470/F |
| Scale / Schaal: 1:5 | Drafted / Getekend | Verified / Nagezien | Confirmed / Bevestigd |  |
| Date / Datum | 29/07/2008 | 21/10/2015 | 21/10/2015 | |
| Init. / Init. | SD | TT | PJ | |