



A	SEAT DEPTH	380, 400, 420, 440, 460, 480, 500
B	SEAT WIDTH	390, 450, 510
C	SEAT ANGLE	0°, 2°, 4°, 6°
D	FOOTREST LENGHT	347, 389, 424, 449, 474,
E	ARMREST HEIGHT	245, 260, 275, 290, 305, 320
F	WIDTH OF ARMREST	411, 471, 531 width adjustment
G	SEAT HEIGHT	440, 490, 540
H	FOOTPLATE HEIGHT	180
I	BACKREST ADJ	0°, 8°, 16°, 24°, 32°



EP00483 / A

DIMENSION TURIOS BASIC				Plan No.
				EP00483/ A
Scale: 1:10	Drafted	Verified	Confirmed	
Date	16/01/2020	12/02/2020	12/02/2020	
Initials	AS	DB	TK	